

# PRACTICE SAFE SPIROMETRY


## Infection Protection For Your Patient, Your Environment ...and You!

SDI Diagnostics, a leading manufacturer in spirometry supplies, has over 20 years experience in PFT filters.


*Pulmo***GUARD IQ**™

### Midmark IQspiro™ spirometer users— protection for you: Save money!

Hollow pneumotachs do nothing to protect you in your test environment. When patients blow, millions of potentially harmful pathogens are forcefully blown into your space.



For your protection, SDI has designed filters specifically for the IQspiro pneumotach.



#### Introducing the PulmoGuard IQ, a member of our new PulmoGuard family

- Feel safer and save money
- Unique, pleated filter media
- Clear plastic for 100% visual inspection
- Recyclable
- Convenient packaging



**sdi**  
DIAGNOSTICS

The Spirometry Source

  
Made In  
America

# PulmoGUARD IQ™

- Lower cost
- 99.99% Bacterial and viral efficiency
- Infection protection for patients and health care professionals
- Exceeds all ATS/ERS standards for spirometry



## PulmoGuard IQ Reorder Info

Part No.	Filter Description	Qty	Packaging
F0601-025	PulmoGuard IQ (25 PulmoGuard IQ filters)	25	25
F0601-100	PulmoGuard IQ (100 PulmoGuard IQ filters)	100	4x25
F0601-300	PulmoGuard IQ (300 PulmoGuard IQ filters)	300	300
<b>Filter Kits</b>			
F0611-025	PulmoGuard IQ (25 PulmoGuard IQ filters, 25 Kushion Klips)	25	25
F0611-100	PulmoGuard IQ (100 PulmoGuard IQ filters, 100 Kushion Klips)	100	4x25
F0611-300	PulmoGuard IQ (300 PulmoGuard IQ filters, 300 Kushion Klips)	300	300 Bulk



Made in America



Your Dealer Is:



CALL 800-728-8230 or 508-588-8222  
 FAX 508-580-8884 EMAIL [sales@quesetmedical.com](mailto:sales@quesetmedical.com)  
 VISIT US ONLINE AT [www.quesetmedical.com](http://www.quesetmedical.com)



DIAGNOSTICS  
 The Spirometry Source  
 Easton, MA  
[www.sdidiagnostics.com](http://www.sdidiagnostics.com)