

# DD-300 Non-display Doppler

# Newman Medical

Newman Medical is proud to offer a full line of Dopplers with the performance and durability you expect to meet your application - obstetrical or vascular. The **DigiDop 300 non-display** Doppler has unique features not found in other Doppler systems such as Digital Optimized Processing (DOP™). Each DigiDop is designed with the sensitivity you need while also designed to be tough and last through years of everyday office use. All our products are **Made in the USA in Colorado**. Each Newman Medical Doppler and probe must pass a series of comprehensive tests to assure quality and the durability of each product.



## Features:

- ✦ **Best Sound** - Digitally Optimized Processing – DOP™ produces up to 6x signal of other Dopplers
- ✦ **Best Warranty** – 5 year \$0 cost warranty – **ALL** Parts. **ALL** Labor.
- ✦ **Best Support & Service** – Knowledgeable Staff & Guaranteed Service
- ✦ Tough and durable polycarbonate and stainless steel case
- ✦ 'Unbreakable' magnetic probe holder
- ✦ All probes are interchangeable
- ✦ Each main unit comes with choice of one probe
  - ✦ 2MHz, 2MHz WP, 3MHz, 3MHz WP, 5MHz, 8MHz
- ✦ User replaceable cord
- ✦ Powered by 3 AA batteries offer extended battery life (900 mins)
- ✦ Dimensions (h x w x l): 150 x 65 x 35 mm (6.0 x 2.5 x 1.25 in)
- ✦ Weight: 340 grams (12 ounces)



 **Queset MEDICAL™**  
P.O. Box 1287  
Brockton, MA 02303

800-728-8230 / 508-588-8222  
Fax 508-580-8884  
[www.quesetmedical.com](http://www.quesetmedical.com)

NM0004B