

**The Leader In  
PFT Filters!**


# PulmoGUARD™

*Fits most PFT devices without adapters*

During a Pulmonary Function Test, bacteria and viruses can be expelled, posing serious risk to patients and health care workers alike. With a cross-contamination efficiency of 99.999%, Pulmoguard is safe, highly efficient...and very affordable.

## PulmoGUARD™

INFECTION PROTECTION



- Meets ATS air resistance standards
- Features dead space of 60cc
- Pleated media for better functionality
- Compact, non-intimidating design
- Transparent for easy, 100% visual inspection
- Fits most spirometers with no adapter

*Please see back for Technical Specifications and uses.*

### PulmoGuard Part Numbers & Packaging

Part No.	Description	Qty	Packaging
F0100-050	PulmoGuard Filter	50	Box
F0100-100	PulmoGuard Filter	100	Box
F0100-200	PulmoGuard Filter	200	4 Boxes/50
F0100-300	PulmoGuard Filter	300	Box

*Please see back for Kit options*

## PulmoGUARD II™

INFECTION PROTECTION



All the benefits of the original Pulmoguard **PLUS –**

- Fits many PFT instruments with no adapter
- Features dead space of 64cc

*Please see back for Technical Specifications and uses.*

### PulmoGuard II Part Numbers & Packaging

Part No.	Description	Qty	Packaging
F0200-050	PulmoGuard II Filter	50	Box
F0200-100	PulmoGuard II Filter	100	Box
F0200-200	PulmoGuard II Filter	200	4 Boxes/50
F0200-300	PulmoGuard II Filter	300	Box

*Please see back for Kit options*



CALL **800-728-8230** or 508-588-8222  
 FAX 508-580-8884 EMAIL [sales@quesetmedical.com](mailto:sales@quesetmedical.com)  
 VISIT US ONLINE AT [www.quesetmedical.com](http://www.quesetmedical.com)

## PULMOGUARD TECHNICAL SPECIFICATIONS

PulmoGuard fits the following instruments without an adapter:

*Adapters available for instruments not listed.*

**Note:** Use PulmoGuard II for all SDI models.

**Assess** - Peak Flow Meter

**Burdick** - Spirosense

**Circadian** - Models (if using black Tamarac pneumotach, need 7975)

**Collins** - Survey, Survey II, Survey III, Survey PLUS, Eagle, Eagle II

**Cybermedics** - Spinnaker TL, SL, CM3, CM310 (if white or black cone use 7980)

**Cybermedics** - Moose

**Godart** - All Models

**Gould** - 5004

**Healthscan** - Assess Peak Flow Meter

**Jaeger** - All Models except Spiro Pro

**Jones** - All Models except Satellite (small end of filter)

**Koko** - All Models

**Med Graphics** - 1070, 1085 (w/Prevent, use PG II; if no Prevent, 7970-spirometry, 7961-vol/diff)

**Med Science** - 570, 6257

**Med Science** - 3000, 4000 (spirometry-no adapt./vol-7961/diff-7960B)

**Med Systems** - Spirolite 101 & 201 (small end of filter), 343

**New Med** - 7000 I (Spirobyte)

**PK Morgan** - Benchmark (spirometry-7978, vol&diff-7960B), Flexi Flow (small end of filter), Trans Flow

**Pneumedics** - All Models (small end of filter)

**Riko** - AS300

**S&M** - CF-150B, GS520, Keystone III, Portascreen, Pulmoscreen

**ViaSys/Sensormedics** - Horizon, Ohio, 922, 2130, 2200, 2400, 2450, 2800, 6200, 2100 (spirometry-7961, vol-7960B, diff-7960B), Vmax (can use PGII w/7976 or 7977)

**Spirolink** - All Models

**Spirometrics** - LTE, 2500 LTE

**Spiroscreen** - All Models (small end of filter)

**Spirotech** - All Models

**Tamarac** - All Models except Presto Flash

## PULMOGUARD II TECHNICAL SPECIFICATIONS

PulmoGuard II fits the following instruments without an adapter:

*Adapters available for instruments not listed.*

**Breon** - All Models

**CDX** - Spiro 110 (remove black fitting), Spiro 850

**Cavitron/Donti** - All Models

**Clement-Clarke** - Mini-Wright Peak Flow Meter

**Collins** - All DS series pulmonary labs with UBV (use 7982 adapter for hose), CPL, GS, SurveyTach, Eaglet, Hawk

**Cosmed** - All Models

**Creative Biomedics** - All Models

**Futuremed** - Discovery, Spirobank, Spirovision III

**Medical Equip Designs** - All Models (remove black fitting)

**Med Graphics** - preVent, PF/DX (needs to be reprogrammed if no filter before)

**Micro Medical/Direct** - All Models

**MIR** - All Models

**Multispiro** - All Models except PC (remove black fitting)

**nSpiro** - Wright Peak Flow Meter

**Ohio** - All Models

**PK Morgan** - 232

**Puritan-Bennett** - PS600, 900, 900L, VS400

**Riko** - All Models except AS300

**S&M** - Keystone

**SDI** - SBOffice, SBG, spirotel, Spirolab, Spirolab II, Astra series with 7983 adapter

**Smiths Medical** - All models

**Spirometrics** - All Models except LTE, 2500 LTE

**VacuMed** - Discovery, Micro, MicroPlus, MicroGP/DL, MicroLoop, Pony, Spiroscreen, Wizard

**ViaSys** - MicroLoop, MicroLab

**Vitalograph** - All Models Wright Professional Peak Flow Meter except Presto Flash

## PulmoGuard Kits

Part No.	Description	Qty	Packaging
F1010-050	PulmoGuard Filter + Cardboard Mouthpiece	50	Box
F0101-100	PulmoGuard Filter + Cardboard Mouthpiece	100	Box
F0102-040	PulmoGuard Filter + Cardboard Mouthpiece + The Klip	40	Box
F0102-080	PulmoGuard Filter + Cardboard Mouthpiece + The Klip	80	Box
F0103-050	PulmoGuard Filter + The Klip	50	Box
F0103-100	PulmoGuard Filter + The Klip	100	Box
F0104-040	PulmoGuard Filter + Comfit-D Disposable Mouthpiece	40	Box
F0104-080	PulmoGuard Filter + Comfit-D Disposable Mouthpiece	80	Box
F0105-040	PulmoGuard Filter + Comfit-D Disposable Mouthpiece + The Klip	40	Box
F0105-080	PulmoGuard Filter + Comfit-D Disposable Mouthpiece + The Klip	80	Box

## PulmoGuard II Kits

Part No.	Description	Qty	Packaging
F0201-050	PulmoGuard II Filter + 1" O.D. Mouthpiece	50	Box
F0201-100	PulmoGuard II Filter + 1" O.D. Mouthpiece	100	Box
F0204-050	PulmoGuard II Filter + 1-1/16" O.D. Mouthpiece	50	Box
F0204-100	PulmoGuard II Filter + 1-1/16" O.D. Mouthpiece	100	Box
F0202-040	PulmoGuard II Filter + 1" O.D. Cardboard Mouthpiece + The Klip	40	Box
F0202-080	PulmoGuard II Filter + 1" O.D. Cardboard Mouthpiece + The Klip	80	Box
F0203-050	PulmoGuard II Filter + The Klip	50	Box
F0203-100	PulmoGuard II Filter + The Klip	100	Box
F0205-040	PulmoGuard II Filter + 1-1/16" O.D. Cardboard Mouthpiece + The Klip	40	Box
F0205-080	PulmoGuard II Filter + 1-1/16" O.D. Cardboard Mouthpiece + The Klip	80	Box
F0208-040	PulmoGuard II Filter + Comfit-D Disposable Mouthpiece	40	Box
F0208-080	PulmoGuard II Filter + Comfit-D Disposable Mouthpiece	80	Box
F0209-040	PulmoGuard Filter + Comfit-D Disposable Mouthpiece + The Klip	40	Box
F0209-080	PulmoGuard Filter + Comfit-D Disposable Mouthpiece + The Klip	80	Box